

TECHNICAL DATA SHEET

Part#: 159B-RTWX

Tweezers - Replaceable Acetal Tips for 159B-RTW

The 159B-RTWX replacement tips for 159B-RTW contain two tips and 4 screws. Acetal has outstanding resistance to chemicals and solvents. These Acetal tips are electrically insulative and resistant to temperatures up to 185°F. The bodies are made of anti-magnetic Neverust® stainless steel (AISI 316L). Excellent for cleanroom applications and for avoiding the marring or scratching of items being handled.

Precision Level
 Material
 Surface Resistivity
 Shape
 ★★★
Acetal
Insulative
Straight

• **Group** Replaceable Tip Tweezers

• Weight (OZ) 0.5

• Tip Size (CTM) .04" x .038" (1mm x .95mm)

• Tip Angle Straight

Tip LengthBody Width

Body Height One Leg n/a
 OAL n/a
 Country of Origin ltaly

UPC
 Schedule B
 Harmonization Code
 ECCN
 UNSPSC Code
 662847004537
 8203.20.20.00
 8203.20.20
 EAR99
 27112105





