

TECHNICAL DATA SHEET

Part#: 181-10

Scalpel Blade #10 Sterile for 181-SE Handle - Round -★★- SS

The 181-10 two star sterile scalpel blade is packaged in individual foil packs. The #10 blade has a curved cutting edge. This blade works with a #3 scalpel handle.

• Precision Level	★★
• Material	Carbon Steel
• Surface Resistivity	Conductive
• Shape	Round
• Group	Cutting Tools
• Blade Angle	
• Blade Length	
• OAL	n/a
• Weight (OZ)	0.1
• UPC	
• Schedule B	8211.95.50.00
• ECCN	EAR99
• UNSPSC Code	
• Country of Origin	Pakistan