

TECHNICAL DATA SHEET

Part#: 181-22

Scalpel Blade Sterile for 181A-SE Handle - Round -★★- SS

The 181-22 two star sterile scalpel blade is packaged in individual foil packs. The blade is similar to the #10 blade but has a curved cutting edge and a flat unsharpeed back edge. This blade works with a 181A-SE scalpel handle.

Precision Level ★★

Material Carbon SteelSurface Resistivity Conductive

Shape RoundGroup Cutting Tools

Blade AngleBlade Length

• OAL n/a • Weight (OZ) 0.1

• UPC 662847012570 • Schedule B 8211.95.50.00

• ECCN EAR99

• UNSPSC Code

• Country of Origin Pakistan



