

TECHNICAL DATA SHEET

Part#: 181A-SE

Scalpel Handle for 181-22 - Straight -★★- SS

The 181A-SE two star scalpel handle works with the 181-22 blades.

• Precision Level	★★
• Material	SS (AISI 440)
• Surface Resistivity	Conductive
• Shape	Straight
• Group	Cutting Tools
• Blade Angle	n/a
• Blade Length	
• OAL	5.5" (137.5mm)
• Weight (OZ)	0.85
• UPC	662847012617
• Schedule B	8211.95.50.00
• ECCN	EAR99
• UNSPSC Code	
• Country of Origin	Pakistan

