

## TECHNICAL DATA SHEET

Part#: 24-4-SE

Forceps - Laboratory - Curved -- SS - Serrated Tips/Grips

*The 24-4-SE eye dressing forceps are thumb forceps used for holding dressing materials during ophthalmic procedures. This product is straight with extra delicate, serrated tips and a length of 4 inches. The grips are serrated for positive control.*

- Precision Level ★
- Material SS (AISI 440)
- Surface Resistivity Conductive
- Shape Curved
- Group Tweezers (Forceps)
- Weight (OZ) 1
- Tip Size (CTM) .02" x .030" (.5mm x .75mm)
- Tip Angle 45°
- Tip Length
- Body Width .276"(6.9mm)
- Body Height One Leg .048"(1.2mm)
- OAL 4"(100mm)
- Country of Origin Pakistan
- UPC 662847011962
- Schedule B 8203.20.20.00
- Harmonization Code 8203.20.20
- ECCN EAR99
- UNSPSC Code 27112105

