

TECHNICAL DATA SHEET

Part#: 272

Scissors - Straight -★★- SS - Blade Length 1.5"

The 272 is a dressing scissors featuring one blunt tipped blade and one pointed blade. This tool is made of surgical stainless steel and is 4.5" long with cutting blades 1.5" long.

Precision Level

**

Material

SS (AISI 440)

Surface Resistivity

Conductive Straight

ShapeGroup

Cutting Tools

Blade Angle

1.5"(37.5mm)

Blade Length

4.56"(114mm)

OALWeight (OZ)

0.72

• UPC

662847011856

Schedule B

8203.00.90.01

• ECCN

EAR100

• UNSPSC Code

44121618

Country of Origin

Pakistan



