

TECHNICAL DATA SHEET

Part#: 275

Scissors - Straight -★★- SS - Blade Length 1.1"

**

SS (AISI 440)

Conductive

Cutting Tools

Straight

The 275 'Iris' scissors is designed for very fine dissection work. This tool is made of surgical stainless steel and is 4" long with cutting blades 1.1" long.

Precision Level

Material

Surface Resistivity

• Shape

• Group

• Blade Angle

• Blade Length 1.1"(27.5mm)

• OAL 4.5" (112.5mm)

• Weight (OZ) 0.629

• **UPC** 662847011887

• Schedule B 8203.00.90.00

• ECCN EAR99

• UNSPSC Code 44121618

• Country of Origin Pakistan





