

## TECHNICAL DATA SHEET

Part#: 275A

---

### Scissors - Curved -★★- SS - Blade Length 1"

The 275A 'Iris' scissors is designed for very fine dissection work. It features a 10° curved blade for visibility. This tool is made of surgical stainless steel and is 5" long with cutting blades 1" long.

• Precision Level	★★
• Material	SS (AISI 440)
• Surface Resistivity	Conductive
• Shape	Curved
• Group	Cutting Tools
• Blade Angle	10°
• Blade Length	1" (25mm)
• OAL	5" (125mm)
• Weight (OZ)	0.555
• UPC	662847011870
• Schedule B	8203.00.90.00
• ECCN	EAR99
• UNSPSC Code	44121618
• Country of Origin	Pakistan

