

TECHNICAL DATA SHEET

Part#: 31-SE

Clamp - Mosquito - Curved -- SS

The 31-SE 30° curved mosquito clamp is used to clamp small blood vessels. The locking mechanism allows the clamp to hold the vessel securely. The serrated jaws allow for positive gripping.

Precision Level ★★

Material SS (AISI 440)
Surface Resistivity Conductive
Shape Curved

• **Group** Hemostats/Clamps

Weight (OZ)

• Tip Size (CTM)

• Tip Angle 30°

Tip Length

Body Width 3.77"(94.25mm)
Body Height One Leg
OAL 5.125"(128mm)
Country of Origin
UPC 62847011849
Schedule B 9018.90.80
Harmonization Code 9018.90.80

ECCN EAR99UNSPSC Code 42151685



