

TECHNICAL DATA SHEET

Part#: 355

Scissors - Medical Grade - Sharp Point - ★★★ - SS - Blade Length .885"

The 355 extremely precise tempered stainless steel (AISI 420) scissors are ideal for delicate cutting applications. They provide a comfortable grip with a contemporary design. The sleek, twisted frame gives them a professional and unique look. They are acid proof, rust resistant and autoclavable.

Precision Level

Material

SS (AISI 420)

Surface Resistivity

Conductive

Shape

Angled

• Group

Cutting Tools

• Blade Angle

.885" (22mm)

Blade Length

3.933" (98.33mm)

OALWeight (OZ)

0.56

• UPC

662847013515

Schedule B

8213.00.90.00

• ECCN

EAR99

• UNSPSC Code

44121618

Country of Origin

Italy

