

TECHNICAL DATA SHEET

Part#: 356

Scissors - Medical Grade - Sharp Point - ★★★ - SS - Blade Length .865"

The 356 extremely precise tempered stainless steel (AISI 420) scissors are ideal for delicate cutting applications. They provide a comfortable grip with a contemporary design. The sleek, twisted frame gives them a professional and unique look. They are acid proof, rust resistant and autoclavable.

Precision Level

Material

SS (AISI 420)

Surface Resistivity

Conductive Angled

• Shape

Cutting Tools

Group

Cutting 100is

• Blade Angle

.865" (22mm)

Blade LengthOAL

4.027" (101mm)

• Weight (OZ)

0.7

• UPC

662847013522

Schedule B

8213.00.90.00

ECCNUNSPSC Code

EAR99

• Country of Origin

44121618

Italy



