

TECHNICAL DATA SHEET

Part#: 357

Scissors - Medical Grade - Sharp Point - ★★★ - SS - Blade Length 1.11"

The 357 extremely precise tempered stainless steel (AISI 420) scissors are ideal for delicate cutting applications. They provide a comfortable grip with a contemporary design. The sleek, twisted frame gives them a professional and unique look. They are acid proof, rust resistant and autoclavable.

Precision Level

. .

Surface Resistivity

ShapeGroup

Material

• Blade Angle

• Blade Length

• OAL

• Weight (OZ)

• UPC

• Schedule B

• ECCN

• UNSPSC Code

Country of Origin

SS (AISI 420)

Conductive

Straight

Cutting Tools

1.1125" (28mm)

3.836" (95.90mm)

0.52

662847013539

8213.00.90.00

EAR99

44121618

Italy





