

TECHNICAL DATA SHEET

Part#: 358

Scissors - Medical Grade - Medium Sharp Point - ★★★ - SS - Blade Length 1.786"

The 358 extremely precise large relieved blade tempered stainless steel (AISI 420) scissors are ideal for delicate applications. Their sharp tips are designed for a wide variety of applications, while providing a superior cut. The sleek, twisted frame gives them a professional and unique look. They are acid proof, rust resistant and autoclavable.

Precision Level

Material

Surface Resistivity

Shape

• Group

• Blade Angle

• Blade Length

• OAL

• Weight (OZ)

• UPC

Schedule B

• ECCN

• UNSPSC Code

Country of Origin

SS (AISI 420)

Conductive

Straight

Cutting Tools

1.786" (44.65mm)

5.375" (134.38mm)

1.57

662847013546

8213.00.90.00

EAR99

44121618

Italy

