

TECHNICAL DATA SHEET

Part#: 36A-SE

Hemostats - 5.5" 25℃urved - SS

The 36A-SE stainless steel clamp is used to clamp small blood vessels. The locking mechanism allows the clam to hold the vessel securely. It has serrated jaws for positive gripping.

Precision Level ★★

Material SS (AISI 440)
Surface Resistivity Conductive
Shape Curved

• Group Hemostats/Clamps

• Weight (OZ) 1.12

• Tip Size (CTM)

• Tip Angle 25°

Tip Length

Body Width
Body Height One Leg
OAL
Country of Origin
UPC
Schedule B
0.26" (6.5mm)
0.172"(4.3mm)
5.5" (137.5mm)
Pakistan
662847012679
9018.90.80.01

• Harmonization Code 9018.90.80 • ECCN EAR99 • UNSPSC Code 42151685





