

## **TECHNICAL DATA SHEET**

Part#: 40-SE

Hemostats - 8" 30℃urved - SS

The 40-SE stainless steel locking hemostat with 30° curved serrated jaws is clean room safe. The overall length is 8".

Precision Level

Material SS (AISI 440)
Surface Resistivity Conductive
Shape Curved

• **Group** Hemostats/Clamps

• Weight (OZ)

• Tip Size (CTM)

• Tip Angle 30°

Tip Length

Body Width .283" (7.08mm)
Body Height One Leg
OAL 8" (200 mm)
Country of Origin
UPC 662847010736
Schedule B 9018.90.80.00

Harmonization Code 9018.90.80
ECCN EAR99
UNSPSC Code 42151685





