

## **TECHNICAL DATA SHEET**

Part#: 5-SA-DC

Tweezers - 4 Star Straight Tapered Ultra Fine Point - Anti-Mag. SS - Diamond Coated

The 5-SA-DC straight tapered ultra-fine pointed tweezers is primarily used for microscopy. The style 5 has the finest points available. The overall length is 4.375". It is made of anti-magnetic Neverust® stainless steel. The diamond coated tips are biocompatible, chemically inert to 350°C, prevent particle shedding, cleanroom safe, and protect the tips from wear while providing elasticity.

Precision Level ★★★★

• Material Anti-Mag. SS (316L)

Surface Resistivity
Shape
Group
Conductive
Straight Tapered
Exotic Metals/Coatings

Weight (OZ)41

• Tip Size (CTM) .004" x .002" (.1mm x .05mm)

• Tip Angle Straight

Tip Length

Body Width .375" (9.38mm)
Body Height One Leg
OAL 4.375" (109mm)
Country of Origin
UPC 662847013133
Schedule B 8203.20.20.00

Schedule B
Harmonization Code
ECCN
UNSPSC Code
8203.20.20
EAR99
27112105



