

TECHNICAL DATA SHEET

Part#: 5A-S-SE

Tweezers - 1 Star Offset Ultra Fine Point - SS

The 5A-S-SE 10° offset angle tapered ultra-fine point tweezer is primarily used for microscopy. The offset angle offers easy access and is designed for good visibility. The overall length is 4.5". It is made of stainless steel.

Precision Level

• Material SS (AISI 440)

Surface Resistivity ConductiveShape Offset


• **Group** Conventional Style Tweezers

Weight (OZ) 1
Tip Size (CTM) N/A
Tip Angle 10°

• Tip Length

Body Width
Body Height One Leg
OAL
Country of Origin
0.393" (9.83mm)
0.06" (1.5mm)
4.5" (112.5mm)
Pakistan

Country of Origin
 UPC
 Schedule B
 Harmonization Code
 ECCN
 UNSPSC Code
 Pakistan
 8203.20.20.00
 8203.20.20.00
 EAR99
 27112105



C) 0.393" D) 4.50" E) 0.060"





5A-S-SE

