

TECHNICAL DATA SHEET

Part#: AF-2-SE

Forceps - Laboratory - Straight .065" x .071" Serrated Tips - SS

This Adson tissue 4.75" medium forceps (1 X 2 teeth) is used for holding and manipulating delicate tissues. These forceps have a wide, flat thumb grasp area that is commonly serrated. The satin finish reduces glare.

- Precision Level ★
- Material SS (AISI 440)
- Surface Resistivity Conductive
- Shape Straight Tapered
- Group Tweezers (Forceps)
- Weight (OZ) 1
- Tip Size (CTM) .065" x .071" (1.63mm x 1.78mm)
- Tip Angle Straight
- Tip Length
- Body Width 0.44"(11mm)
- Body Height One Leg 0.04"(1mm)
- OAL 4.75" (118.75mm)
- Country of Origin Pakistan
- UPC 662847011948
- Schedule B 8203.20.20.00
- Harmonization Code 8203.20.20
- ECCN EAR99
- UNSPSC Code 27112105

