

TECHNICAL DATA SHEET

Part#: F-2-SE

Forceps - Laboratory - Curved .02" x .03" Serrated Tips - SS

The F-2-SE eye dressing forceps are thumb forceps used for holding dressing materials during ophthalmic procedures. This product is straight with extra delicate, serrated tips and a length of 4 inches. The serrated grips facilitate better control of the tool.

Precision Level

Material SS (AISI 440)
Surface Resistivity Conductive
Shape Curved

• **Group** Tweezers (Forceps)

• Weight (OZ) 0.39

• Tip Size (CTM) .02" x .03" (.5mm x .75mm)

EAR99

• Tip Angle 90°

Tip Length

• ECCN

Body Width 0.306"(7.65mm)
Body Height One Leg
OAL 4"(100mm)
Country of Origin Pakistan
UPC 62847012297
Schedule B 8203.20.20.00
Harmonization Code 8203.20.20

• UNSPSC Code 27112105



