

## **TECHNICAL DATA SHEET**

Part#: F-4-SE

Forceps - Laboratory - Curved .015" x .025" Serrated Tips - SS

The F-4-SE is a 45° half curved 3.75" dressing forceps used for compressing or grasping tissues, handling sterile dressings, and other purposes. The serrated tips give more gripping power, while the half curved tips allow for better visibility.

Precision Level

Material SS (AISI 440)
 Surface Resistivity Conductive
 Shape Curved Tapered
 Group Tweezers (Forceps)

• Weight (OZ) 0.39

• Tip Size (CTM) .015" x .025" (.38mm x .63mm)

• Tip Angle 45°

Tip Length

Body Width 0.411"(10.28mm)
 Body Height One Leg 0.063"(1.58mm)
 OAL 3.75"(93.75mm)
 Country of Origin Pakistan
 UPC 662847012280
 Schedule B 8203.20.20.00

Harmonization Code
ECCN
UNSPSC Code
EAR99
27112105



